





CARSOFT 4	UNIMAT
	
ROTAMAT	BARIA FORM
	

Technical Manual

CARSOFT 4 mattress

UNIMAT mattress

ROTAMAT mattress

BARIA FORM mattress

HALFA mattress

SOFFIA mattress

TABLE OF CONTENTS

TECHNICAL CHARACTERISTICS OF THE DEVICE	- 3 -
COMPATIBLES PROTECTIVE ACCESSORIES	- 5 -
PRECAUTIONS FOR USE	- 5 -
PREVENTIVE MAINTENANCE	- 6 -
CARE AND DISINFECTION	- 6 -
Cleaning the cover	- 6 -
Caution	- 7 -
DISPOSAL	- 7 -

TECHNICAL CHARACTERISTICS OF THE DEVICETechnical characteristics:

Technical characteristic of UNIMAT	
Top layer of high-elasticity foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 1.6 Kpa
Bottom layer of high resilience foam	Polyurethane 37 kg/m ³ – Load-bearing capacity: 3.5 Kpa Fire classification: M4
Insert of high resilience foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 6 Kpa Fire classification: M4
Integral cover in CARTEX	65% Polyurethane – 35% Polyester Weight: 140 g/cm ² Fire classification: crib 5
Integral cover in bi-elastic PHARMATEX	50% Polyurethane – 50% Polyester Grammage: 205 g/cm ² Fire classification: crib 5
Fire resistance test	EN 597-1 & EN 597-2

Technical characteristic of CARSOFT 4	
Top layer of viscoelastic foam	Polyurethane 50 kg/m ³ – Load-bearing capacity: 1.8 Kpa Fire resistance test: EN 1021-1
Bottom layer of high resilience foam	Polyurethane 30 kg/m ³ – Load-bearing capacity: 6.3 Kpa
Insert of high elasticity foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 1.6 Kpa
Integral cover in CARTEX	65% Polyurethane – 35% Polyester Weight: 140 g/cm ² Fire classification: crib 5
Fire resistance test	EN 597-1 & EN 597-2

Technical characteristic of ROTAMAT	
Outer layers of high elasticity foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 1.6 Kpa
Insert of high resilience foam	Polyurethane 37 kg/m ³ – Load-bearing capacity: 3.5 Kpa Fire classification: M4
Integral cover in CARTEX	65% Polyurethane - 35% Polyester Weight: 140 g/cm ² Fire classification: crib 5
Integral cover in bi-elastic PHARMATEX	50% Polyurethane – 50% Polyester Grammage: 205 g/cm ² Fire classification: crib 5
Fire resistance test	EN 597-1 & EN 597-2

Technical characteristic of BARIA FORM

Top layer of viscoelastic foam	Polyurethane 50 kg/m ³ – Load-bearing capacity: 2.5 Kpa Fire classification: M4
Bottom layer of high resilience foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 6 Kpa Fire classification: M4
Middle layer of high resilience foam	Polyurethane 37 kg/m ³ – Load-bearing capacity: 3.5 Kpa Fire classification: M4
Integral cover in CARTEX	65% Polyurethane – 35% Polyester Weight: 140 g/cm ² Fire classification: CRIB 5
Anti-slip undercover	80% PVC – 20% Polyester
3D mesh	50% Polyurethane – 50% Polyester Weight: 205 g/cm ²
Fire resistance test	EN 597-1 & EN 597-2

Technical characteristic of SOFFIA

Top layer of viscoelastic foam	Polyurethane 50 kg/m ³ – Load-bearing capacity: 2.5 Kpa Fire classification: M4
Bottom layer of high resilience foam	Polyurethane 37 kg/m ³ – Load-bearing capacity: 3.5 Kpa Fire classification: M4
Middle layer of high resilience foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 1.6 Kpa
Insert of high resilience foam	Polyurethane 40 kg/m ³ – Load-bearing capacity: 6 Kpa Fire classification: M4
Integral cover in CARTEX	65% Polyurethane – 35% Polyester Weight: 140 g/cm ² Fire classification: CRIB 5
Anti-slip undercover	80% vinyl of polychloride (PVC) – 20% Polyester
Fire resistance test	EN 597-1 & EN 597-2

Technical characteristic of HALFA

Top layer of viscoelastic foam	Polyurethane 50 kg/m ³ – Load-bearing capacity: 2.5 Kpa Fire classification: M4
Bottom layer of high resilience foam	Polyurethane 28 kg/m ³ – Load-bearing capacity: 3.4 Kpa Fire classification: M4
Integral cover in PROMUST	58% Polyurethane – 42% Polyester Weight: 120 g/cm ² Fire classification: CRIB 5
Fire resistance test	EN 597-1 & EN 597-2

COMPATIBLES PROTECTIVE ACCESSORIES

REFERENCES	DIMENSIONAL CHARACTERISTICS
CARSOFT 4 cover	
639212 for CARSOFT 4 mattress – BED 90	88 x 198 x 16 cm
ROTAMAT covers	
04.416DH, 04.416DDH, 04.416TH, 04.416TDH for ROTAMAT mattress - BED 90	88 x 198 x 16 cm
BARIA FORM covers	
04.6090H for BARIA FORM mattress – BED 90	90 x 195 x 20 cm
04.6100H for BARIA FORM mattress – BED 100	100 x 195 x 20 cm
04.6110H for BARIA FORM mattress – BED 110	110 x 195 x 20 cm
04.6120H for BARIA FORM mattress – BED 120	120 x 195 x 20 cm
04.6140H for BARIA FORM mattress – BED 140	140 x 195 x 20 cm
04.6160H for BARIA FORM mattress – BED 160	160 x 195 x 20 cm
UNIMAT covers	
04.15DH, 04.15DDH, 04.15TH, 04.515TDH for UNIMAT mattress – BED 90	88 x 198 x 16 cm
SOFFIA covers	
17.HMS602H for SOFFIA mattress – BED 90	198 x 88 x 16 cm
17.HMS602-10H for SOFFIA mattress – BED 100	198 x 100 x 16 cm
HALFA covers	
17.HMS601H for HALFA mattress – BED 90	196 x 86 x 14 cm
17.HMS601-10H for HALFA mattress – BED 100	196 x 100 x 14 cm



Only protective accessories distributed by WINNCARE France guarantee compatibility with our products. Incompatible accessories may create hazards.

PRECAUTIONS FOR USE

- A mattress should be installed on a flat and high-quality bed frame,
- For mattresses delivered compressed, check the "use by" date on the box. Respect the decompression time of 2-10 minutes before placing the patient on the mattress.
- The use of a support to prevent pressure ulcers/injuries is accompanied by therapeutic education of the person equipped with this device and their family. This must be provided by a health professional working at home or in an institution. A follow-up booklet of recommendations and preventive actions can be created to record the preventive actions.
- It is important to keep excess thickness between the body and the support as low as possible, except for the outer garment and a possible full body diaper. Prefer cotton underwear that is not too tight and, if possible, without seams in the support areas. Do

not use folded towels or sheets, extra pillows, etc. Ensure that there are no foreign objects such as tubes, crumbs, fat, etc.

- For reasons of hygiene, each pressure ulcer/injury prevention aid should be used by one person only.
- Check that the dimensions of the mattress, bed and bed rails match the user's body shape.
- Check that there is no gap between the bottom rail of the bed rail and the top of the mattress. A minimum distance of 22 cm between the top of the bed rail and the top of the mattress when the patient is present is recommended.
- Do not place the mattress on top of another mattress, but on top of a bed base. Check that the base is in good condition, that it is not sagging in the middle.
- Mark the direction of installation of the mattress. The protective strip should be positioned so that it covers the zip, which is not waterproof.
- Write the date of use of the mattress directly on the cover label in the indicated place with a marker.
- Do not cover the mattress with anything other than the cover or protection recommended by the manufacturer. Use a replacement cover recommended by the manufacturer to maintain the characteristics of the mattress.
- Place a cotton sheet over the cover. Do not tuck it in to avoid a hammock effect.
- In case of incontinence, clean the cover and possibly the mattress. When washing, it is advisable to use a spare cover to avoid direct contact with the mattress surface.
- Avoid using sharp objects on the mattress.
- Reinforced surveillance in subjects at risk of falling (patient profile: invalid, agitated, not lucid).

PREVENTIVE MAINTENANCE

Regularly check the cover for holes or signs of erosion.




Contact your dealer in case of abnormal deterioration.



CARE AND DISINFECTION

CLEANING THE COVER

Surface cleaning and disinfection (bioneting). Compatible with steam processes. Machine washable at a recommended temperature of 60°C with a cold-down principle. Open the zip.

The following conditions must be respected for CARTEX, PROMUST and PHARMATEX covers:

	Moderate wash up to 90°C.
	Maximum allowable chlorine concentration of 5000 ppm.
	Tumble dry on low heat.

	Do not iron.
	Do not dry clean.

CAUTION

Do not use corrosive cleaning products or industrial degreasers to avoid damage to the textile.

DISPOSAL

Do not burn or throw away the used product.

Take it to a waste disposal site, which will take care of it.



WINNCARE France

4, Le Pas du Château
85670 SAINT-PAUL-MONT-PENIT
FRANCE

Tél. : +33 (0)4.66.02.15.15
Fax. : 33 (0) 4.66.02.15.00
Email: contact@winncare.fr
www.winncare.com



WINNCARE France

Zone Actipole, 2 Rue de Saint-Coulban
35590 MINIAC-MORVAN – France
FRANCE

