

ACCESSORIES FOR DYNAMIC MATTRESSES

QUICK INFLATION AND DEFLATION PUMP



QUICK INFLATION PUMP VKIT-AXT

Designed for short-term use*, the quick inflation and deflation pump allows quick and easy installation of the AXTAIR dynamic mattresses.

*Designed to inflate 1 or 2 mattresses successively



QUICK INFLATION PUMP VAXT/PGR

Designed for long-term use**, the quick inflation and deflation pump allows quick and easy installation of the AXTAIR dynamic mattresses.

**Designed to inflate many mattresses successively

DESIGNED TO FACILITATE THE SET UP AND REMOVAL OF AXTAIR DYNAMIC MATTRESSES, THESE PUMPS QUICKLY AND SAFELY INFLATE OR DEFLATE MATTRESSES.

- ▶ One pump for inflating and deflating AXTAIR dynamic mattresses.
- ▶ The pump allows inflating or deflating AXTAIR dynamic mattress in less than 3 minutes, saving time when installing or removing the mattress.

DESCRIPTION	ITEM REFERENCE	DIMENSIONS	WEIGHT POMPE	KEY FEATURES
QUICK INFLATION PUMP	VKIT-AXT	35 x 14 x 12 cm	1.54 kg	Inflation/deflation time: < 3 min Inflation pressure: > 125 mBar Sound level (dBA at 1m): 91 dBA Electric classification: Class 1 Electric frequency: 60 Hz Electric tension: 220 V Electric cord: PVC 3.4 M EC certification
QUICK INFLATION PUMP	VAXT/PGR	30 x 10,5 x 10,5 cm	1.23 kg	Inflation/deflation time: < 3 min Inflation pressure: > 56 mBar Sound level (dBA at 1m): 76 dBA Electric classification: Class 2 Electric frequency: 60 Hz Electric tension: 220 V Electric cord: 3 M rubber EC certification

Cleaning instructions

2 years warranty

CABLE PROTECTION



CABLE PROTECTION

- ▶ Allows the attachment of the pump power supply cable along the mattress.
- ▶ The absence of a cable on the ground reduces the risk of falls of caregivers of patients.
- ▶ Keeping the cable in the protection reduces the risk of pinching or shearing.
- ▶ Cable can be easily positioned and removed by means of snap fasteners.

DESCRIPTION	ITEM REFERENCE	KEY FEATURES
CABLE PROTECTION	VAXT/PC	Polychlorure de Vinyle

Cleaning instructions: Use of a surface detergent-disinfectant is approved
Moderate wash up to 90°C
Use of chlorine solution is permitted. 2 years warranty

RED CARRYING BAG



RED CARRYING BAG

- ▶ Facilitates safe removal and transport of soiled dynamic mattresses
- ▶ The red color prevents any mistakes during handling and prevents mixing a clean product (blue bag) with a soiled product (red bag).
- ▶ The red carrying bag can hold a mattress with its pump and its dynamic cushion.

DESCRIPTION	REFERENCE	KEY FEATURES
RED CARRYING BAG	VAXT/SR	Nylon enduction polyuréthane

Cleaning instructions: Moderate wash up to 60°C
Can be tumble-dried at low temperature (air evacuated at 60°C max.) 2 years warranty

FP000049/01_2018

Photos not contractually binding

