



UNIVERSAL HAMMOCK with head support

Unique sling with entire body support and good comfort.

USAGE

- For less active users, who are in seated position for longer periods.
- For transfers from a lying position.
- For transfers in mobile and stationary hoists.

FEATURES

- The sling is transferring user from lying position or from the bed or bath.
- The sling provides good support of the body, which gives the user a feeling of safety during transfers.
- The sling ensures an upright sitting position that contributes to a correct positioning of the user.
- The sling with head support is dedicated for users, who need extra support around the upper body and head.
- The sling is available with commode opening, which eases the carer's tasks.
- The sling is made of polyester mesh, which is perfect for bathing.
- Universal Hammock has unsplit leg straps

SPECIFICATIONS

SWL: 275 kg

Material: Polyester Mesh

Størrelser: XS-S-M-L-XL-XXL

ITEM NUMBERS:

Universal Hammock with head support, polyester mesh:

	Item no.
XS	009-09683
S	009-09684
M	009-09685
L	009-09686
XL	009-09687
XXL	009-09688




SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XS	S	M	L	XL	XXL
Y (hip measurement)	37 cm	42 cm	48 cm	52 cm	68 cm	78 cm
Z (measurement from head to lower back)	75 cm	80 cm	91 cm	100 cm	100 cm	100 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winncare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

FOR MORE INFORMATION VISIT

winncare.com

E-mail: info@winncare.dk

Phone: +45 70273720

WIN  CARE